

Quality controls, and appropriate calibration, allow you to measure with...

CONFIDENCE.

In the cardiac catheterization lab it is critical. Waters' Oxicom® whole blood oximeters provide cardiologists with critical information on the oxygen saturation of a patient's blood. The Oxicom series is quick, accurate, and reliable, with over a half-million saturation samples measured every year.

Quality Control – daily

Simple, daily use of three QC filters to guide you through the quality control process. When prompted, simply insert the appropriate QC filter.

Quality Control – weekly

Weekly use of liquid controls are important for calibration verification of the Oxicom. Liquid controls are available in three levels – for monitoring performance of significant points within the clinical range.



Liquid controls are packaged 9 per box.
(3 of each level per box.)
Order number 0034091.000.

One Box: \$495.00 ea
Two Boxes: \$445.00 ea

High (red) Mid (yellow) Low (blue)

On-line Result Registration

Visit watersmed.com and click on the link in the lower left corner of the homepage to register your results.